

Hematocrit centrifuge

TG-12



Order No.:010101200

Features

- >TG-12 is for testing hematocrit
- >Brushless motor. Rapid acceleration and deceleration
- >Electronic lid lock. Cannot open lid when rotor is spinning
- >Microprocessor control. High quality steel housing. Stainless steel chamber
- >LED display important parameters in real time and real values

TECHNICAL DATA

Model	TG-12
Max Speed	12000rpm
Max RCF	15800xg
Capacity	24pcs capillary tube
Timer range	0-99m59s
Speed accuracy	±20rpm
RCF and PRM convert function	Yes
Min. Speed Change	10 kinds
Noise	≤55dB(A)
Power Supply	AC 110V/220V 50HZ/60HZ
Dimension	355x270x205mm
Weight	13kg

Rotor	Max speed	Max capacity	Max RCF
 Order No.:10002400017	12000rpm	24pcs capillary tube	15800xg