

Gel card centrifuge

TD-4K



Order No.:0101230000

Features

- >LED screen shows RPM/RCF and Time
- >Low noise,without carbon powder pollution
- >Variable frequency motor drive, microprocessor control
- >Electronic lid lock to ensure safety
- >Steel housing;304 stainless steel chamber
- >Professional program design,it can run directly without setting paramters

TECHNICAL DATA

Model	TD-4K
Max Speed	3840rpm
Max RCF	1790xg
Capacity	24/12 Gel card
Noise	≤60dB(A)
Power Consumption	100w
Dimension	375x300x260mm
Weight	17kg

Centrifugal procedure	Max Speed	Max RCF
0-2Min	900rpm	100×g
2-5Min	1500rpm	280×g



24 Gel card

Order No.:0103500200



12 Gel card

Order No.:0103500100