

Ventilated | Low speed | Benchtop

TD-4Z



Order No.:0101010000

The TD-4Z is a practical centrifuge for small laboratories and clinics. Standard fixed angle rotor 12x20ml is included to hold various blood tube and 15ml tube

Features

- >LED screen shows RPM/RCF and Time
- >RCF can be set directly
- >1min to 99min timer range
- >User can reset parameters when rotor is spinning
- >Variable frequency motor drive
- >Electronic lid lock to ensure safety
- >Steel housing;304 stainless steel chamber
- >Autoclavable aluminum alloy fixed angle rotor
- >Compatible with blood tube and 15ml tube
- >Can run 12 tubes at a time



Order No.:0494500100



Order No.:0494500100



Order No.:0494501200

TECHNICAL DATA

Model	TD-4Z
Max Speed	4200rpm
Max RCF	2680xg
Capacity	4x50ml;12x15ml tube/ blood tube
Speed Accuracy	±20rpm
Time Setting	1min ~ 99min
Noise	≤55dB(A)
Power Consumption	200w
Dimension	440x340x290mm
Weight	17kg