

Mini | Portable | Ventilated

Mini-6K/7K/10K



Order No.:0101350000(Mini-6K)
Order No.:0101360000(Mini-7K)
Order No.:0102180000(Mini-10K)

Features

- >Compact design, light weight
- >Adjustable run time,1s-99m59s
- >800rpm to highest rpm is adjustable
- >Automatically stop spinning when lid is open
- >BLDC motor drive, stable and low noise
- >Can stop spin in 3 seconds, saving experiment time
- >Wide range power input:100v to 240v
- >Three-in-one rotor:8x1.5/2ml+6x0.5ml+2x8x0.2ml
no need to change rotor
- >Optional rotor for 6 PCR strips



Order No.:0494010100



Order No.:0494010400

TECHNICAL DATA

Model	Mini-6K/Mini-7K/Mini-10K
Max Speed	6000rpm/7500rpm/10000rpm
Max RCF	2200xg/3450xg/6100xg
Capacity	Three in one rotor:8x1.5/2ml+6x0.5ml+2x8x0.2ml Rotor for 6 PCR strips(Mini-6k,Mini-7k)
Time Setting	1s ~ 99min59s
Noise	≤52dB(A)/≤52dB(A)/≤57dB(A)
Power Supply	AC 100-240V 50/60HZ
Dimension	200x170x130mm
Weight	1.2kg