PRODUCT SPECIFICATION

EDTA K2/K3 Vacuum Blood Collection Tube

PURPLE CAP | JUNNUO | MODEL CY

Compiled from the manufacturer's published product data



ADDITIVE	CAP	MATERIAL	VOLUME RANGE
EDTA K2 / K3	Purple	PET / Glass	1-10 mL*

01 / OVERVIEW

Product overview

A vacuum blood collection tube with a purple stopper and EDTA K2 / EDTA K3 additive, listed by the manufacturer for clinical hematology testing. The tube is described as available in PET or glass and is intended to collect a predetermined blood volume.

02 / PRODUCT DATA

Technical specifications

Specification	Manufacturer-published information
Brand	JUNNUO
Product name	K2/K3 EDTA Vacuum Blood Collection Tubes
Model	CY
Additive	EDTA K2 / EDTA K3 (exact variant to confirm)
Material	PET or glass
Tube dimensions	13 x 75 mm; 13 x 100 mm; 16 x 100 mm
Available volume	1-10 mL overall range; size-specific mapping not itemized*
Country / region of origin	Shandong, China
Device classification	Class I (as stated by manufacturer)

INTENDED APPLICATIONS

Clinical hematology testing, including:

- Routine blood examination
- Hematocrit and hemoglobin
- Blood cell morphology / count
- Sickling test

** The selected product page gives an overall 1-10 mL range but does not map individual draw volumes to each size.*

Source: manufacturer product page. This is a compiled data sheet, not a manufacturer-issued IFU.

PRODUCT SPECIFICATION

Technical details & supply information

JUNNUO | EDTA K2/K3 VACUUM BLOOD COLLECTION TUBE

03 / HANDLING DATA

Operating and storage data

Parameter	Manufacturer-published information
Centrifugation speed	3500-4000 rpm
Centrifugation time	5 min
Recommended storage	4-25 °C
Shelf life	1 year (selected product page)
Processing field	OZONE (manufacturer wording; scope not clarified)



Product image supplied in the reference product sheet.

04 / CONFIGURATIONS

Available tube dimensions

Tube dimensions	Volume information	Note
13 x 75 mm	1-10 mL overall range	Individual mapping not listed
13 x 100 mm	1-10 mL overall range	Individual mapping not listed
16 x 100 mm	1-10 mL overall range	Individual mapping not listed

05 / SUPPLY

Packaging information

INNER PACK

100 pcs / tray

Manufacturer-published

CARTON OPTIONS

1200 or 1800 pcs

Confirm option for selected size

MODEL

CY

JUNNUO

DOCUMENTATION NOTES

Confirm these points against the current product label and manufacturer documentation:

- **Exact additive:** Confirm whether the supplied version contains EDTA K2 or EDTA K3.
- **Draw volume by size:** The selected page does not provide an individual volume-to-dimension mapping.
- **Shelf life and processing:** The selected page lists a 1-year shelf life and “OZONE” under a processing field; confirm both against the current label / IFU.
- **Use parameters:** Verify centrifugation conditions against the instructions for the exact supplied variant and laboratory protocol.

Manufacturer source: Chengwu Medical / JUNNUO product page.

Compiled from publicly available product information; not a manufacturer-issued or certified technical file.

For procurement, request the current technical sheet, instructions for use, and lot-specific documents from the supplier.